





LED X-Ray Film Illuminators

India's slimmest medical X-Ray film illuminator

- ❖ World's slimmest Film Illuminator
- ❖ Elegant and Compact
- ❖ Flicker-free hence no eye-fatigue
- ❖ Individual ON/OFF control
- ❖ Individual brightness control
- ❖ Max brightness : 5000cd/m2
- ❖ Individual film sensor
- ❖ LED light source
- ❖ 10000 lux, adjustable
- ❖ LED backlit panel
- ❖ 60,000 hours LED lamp life
- ❖ 95% uniform light distribution
- ❖ Color temperature : 5600K
- ❖ Power consumption : approx 25W*



Available Sizes:

- ❖ Single film - 14" X 17"
- ❖ Two films - 28" X 17"
- ❖ Three films - 42" X 17"
- ❖ Four films - 56" X 17"

Authorised Business Associate:



SINCE 1995

[]

[]

Avantec Medical Systems (P) Ltd.,

No. 76, 7th Street, Porur Gardens Phase-1, Chennai-600 095.

Tel: +91 81447 30553/54/55, Email: 724365@avanttec.net

www.avanttec.net